

The Eckart Klobe Vacuum Bell for lifting of funnel chest

The alternative to surgery



Case example 2251, male
13 years, 171 cm, 49 kg

Left:
2.2 cm before treatment (horizontally)
3.5 cm before treatment (vertikally)

Right:
After one month

The Eckart Klobe Vacuum Bell for lifting of funnel chest

The alternative to surgery

- * **Not surgical**
- * **Not invasive**
- * **Performed during daily life activities**
- * **Applicable until middle age**
- * **Functions similar to braces used for teeth**



The 5 different Vacuum Bell Types

The Large Vacuum Bell (26 cm). For persons taller than 1.70 m
The Small Vacuum Bell (19 cm). For persons taller than 1.40 m
The Small Vacuum Bell Type Bodybuilder (19 cm). For bodybuilders and special applications. For persons taller than 1.40 m
The Small Vacuum Bell Type Women (19 cm). For females after beginning of the breast development and taller than 1.40 m
The Mini Vacuum Bell (16 cm). For children taller than 1.05 m

The Eckart Klobe Vacuum Bell for lifting of funnel chest

The alternative to surgery



Case example 1157, male
57 years, 190 cm, 80 kg

Left:
4.0 cm before treatment

Right:
After 4 Months depth 2.0 cm

Eckart Klobe
M 7, 9a-10
68161 Mannheim Germany
Phone +49-621-13528
www.trichterbrust.de
webmaster@trichterbrust.de